

Electrodrive patient movers compatibility chart

Electrodrive Patient Movers are designed to be used with hospital beds, stretchers and trolleys.

This compatibility chart lists products that have been tested successfully with one or more of our patient movers and approved for use by Electrodrive.

If you are cannot find the product you are looking for, please contact our Healthcare Manager on 03 9300 8500.



Gzunda (GZ10SL)



Gzunda (GZS)



Powered Fifth Wheel (P5W)

Electrodrive patient movers compatibility chart

COMPATIBILITY CHECK

Manufacturer	Model/name	COMPATIBILITY CHECK			Accessory bracket	Notes
		GZ10SL	GZS	P5W		
ARJO	500	✓	✓	✓		
ARJO	Contoura 480	✓	✓	✓		
ARJO	Contoura 380	✓	✓	✓		
ARJO	Contoura 460	✓	✓	✓		
ARJO	Sorrento	✓	✓	✓		
ARJO	Enterprise 8000	✓	✓	✓		Foot extension
ARJO	Enterprise 5000	✓	✓	✓		Foot extension
ARJO	Contoura 1080	✓	✓	✓		Foot extension
CGM Taiwan	Vanguard Classic	✓	✓	✓	Recommend using with SIE0369 (GZ10SLSG 'Singapore Slimline')	
Engrit	EN71	✓	✓	✓		Hitch from head-end only
Hausted	Unicare 800	X	✓	✓		
Horizon	Electric Powered Stretcher	X	✓	X		Column bed
Hausted	Youth Stretcher	X	✓	X		Column bed
Hendicare	700	✓	✓	✓		
Hill-Rom	Avantgarde 800	X	✓	X		Column bed
Hill-Rom	Avantgarde 1200	X	✓	X		Column bed

Electrodrive patient movers compatibility chart

COMPATIBILITY CHECK

Manufacturer	Model/name	COMPATIBILITY CHECK			Accessory bracket	Notes
		GZ10SL	GZS	P5W		
Hill-Rom	Avantgarde 1600	X	✓	X		Column bed
Hill-Rom	Total Care Connect	X	✓	X		Column bed
Hill-Rom	Medicraft 500 (MC500)	✓	✓	✓		
Hill-Rom	Medicraft 700 (MC700)	✓	✓	✓		
Hill-Rom	Medicraft 200 (MC200)	X	✓	X		
Hill-Rom	FE5000	✓	✓	✓		
Hill-Rom	HLE 8000	✓	✓	✓		
Hill-Rom	Versacare 3200	✓	✓	✓		
Hill-Rom	HR900	✓	✓	✓		Suitable with Gripper
Howard Wright	M7	✓	✓	✓	If using Gripper, must use Howard Wright Gripper SIE0305 (jaw-lock—50 mm size)	Foot extension
Howard Wright	M8	X	✓	X	If using Gripper, must use Howard Wright Gripper SIE0305 (jaw-lock—50 mm size)	Column bed
Howard Wright	M9	✓	✓	✓	If using Gripper, must use Howard Wright Gripper SIE0305 (jaw-lock—50 mm size)	Foot extension/requires 90mm tine
Ikan	Rose Floorline	X	✓	X		Column bed

Electrodrive patient movers compatibility chart

COMPATIBILITY CHECK

Manufacturer	Model/name	COMPATIBILITY CHECK			Accessory bracket	Notes
		GZ10SL	GZS	P5W		
Invacare	1600	✓	✓	✓		
Invacare	Alegio	✓	✓	✓		
Invacare	Aria	X	✓	X		Column bed
Invacare	ScanBed 750	X	✓	X		
Invacare	Etude Basic	✓	✓	✓		
Invacare	Etude Duo	✓	✓	✓		
Invacare	Etude Classic	✓	✓	✓		
Invacare	Octave	✓	✓	✓		
Linet	Eleganza Smart	X	✓	X		Column bed
Linet	Eleganza Standard	X	✓	X		Column bed
Linet	Eleganza XC	X	✓	X		Column bed
Jackson Shand	ECO	✓	✓	✓		
Jackson Shand	Phoenix	✓	✓	✓		
Jomor	LHL272	✓	✓	✓		
Jomor	ULT222	✓	✓	✓		
Jomor	FLO101	X	✓	X		Column bed

Electrodrive patient movers compatibility chart

COMPATIBILITY CHECK

Manufacturer	Model/name	COMPATIBILITY CHECK			Accessory bracket	Notes
		GZ10SL	GZS	P5W		
Joyce	900/E	✓	✓	✓		
Joyce	Centenary	✓	✓	✓		
Joyce	400/E	✓	✓	✓		
Maxi Promotions	Floor Care	X	✓	X		Column bed
Maxi Promotions	Deluxe Regal	✓	✓	✓		
Maxi Promotions	Low Level	✓	✓	✓		
Maxi Promotions	Ultimate Regal	✓	✓	✓		
Prius	Elite	✓	✓	✓		Foot extension
Prius	CM500	✓	✓	✓		
Liftcare	Protean	X	✓	X		Column bed
Steris	Horizon Series	X	✓	X		Column bed
Stryker	S3	X	✓	X		Column bed
Stryker	Eye Surgery Stretcher	X	✓	X		Column bed
Unique Care	Standard High/Low bed	✓	✓	✓		
Volker	S962-2	✓	✓	✓		